

Medisoft 19 System Requirements

NOTE: Medisoft® 19 uses Advantage 11.10. If you are upgrading, you MUST install Advantage 11.10. Once you have installed Advantage 11.10, you will not be able to access an Advantage Server that is prior to Advantage 11. You will be able to access your database as long as your database server is at least Advantage 11.

NOTE: Medisoft 19 is a 32-Bit application, and on a supported 64-Bit platform, the application will run in a 32-Bit mode.

Medisoft 19 Supported Operating Systems

OPERATING SYSTEM	SERVER		WORKSTATION	
	Version of Medisoft 19		Version of Medisoft 19	
	Network Professional	Single-User or Advanced	Client/Server	Single-User or Advanced
Windows Server 2003 SP2 (32 or 64-Bit)	Y	N/A	Y	Y
Windows Server 2008 SP2 or higher (32-Bit)	Y	N/A	Y	Y
Windows Server 2008 R2 (64-Bit)	Y	N/A	Y	Y
Windows Server 2012 (64-Bit)	Y	N/A	Y	Y
Windows 7 Professional/Ultimate (32 or 64-Bit)	Y	Y	Y	Y
Windows 8 Professional/Enterprise (32 or 64-Bit)	Y	Y	Y	Y
Windows Vista Business SP1 or higher (32 or 64 -Bit)	Y	Y	Y	Y
Windows XP Professional SP3 (32-Bit)	Y	Y	Y	Y

Medisoft will not support Windows RT on the ARM processor.

Medisoft Hardware Requirements

Workstation (minimum required)

CPU (Processor)	Intel Pentium 4 2.0GHz or higher
RAM (Memory)	2GB
Storage Space available	10GB* (20GB for Net Pro)
Optical Drive	DVD-ROM (required if installing from a CD)
Network Card (NIC)	100Mbps
Display Monitor	1024x768 (1280x800 for widescreen displays)

Workstation (Recommended)

CPU (Processor)	Intel i3 or higher
RAM (Memory)	4GB
Storage Space available	20GB*
Optical Drive	DVD-ROM (required if installing from a CD)
Network Card (NIC)	1Gbps
Display Monitor	1024x768 (1280x800 for widescreen displays)

Server (Minimum Required)-Network Professional

CPU (Processor)	Intel Core 2 Duo 2.0GHz
RAM (Memory)	4GB
Storage Space available	20GB*
Optical Drive	DVD-ROM (required if installing from a CD)
Network Card (NIC)	100Mbps
Display Monitor	1024x768 (1280x800 for widescreen displays)

Server (Recommended)-Network Professional

CPU (Processor)	Intel i3 or higher
RAM (Memory)	8GB
Storage Space available	40GB*
Optical Drive	DVD-ROM (required if installing from a CD)
Network Card (NIC)	1Gbps
Display Monitor	1024x768 (1280x800 for widescreen displays)

IMPORTANT: You must also have an Internet connection to use Network Professional.

*The amount of storage space needed will increase as your practice database gets bigger..

Server Operating System – Minimum Required

- Windows XP SP3 (32 Bit)
- Windows Server 2003 (Enterprise or Standard editions with latest service packs)

Server Operating System – Recommended

- Windows 7 Professional/Ultimate (32 or 64 Bit)
- Windows Server 2008 R2 (64 Bit)

Server Operating System – Supported

- Windows XP SP3 (32 Bit)
- Windows Vista Business SP1 or higher (32 or 64 Bit)
- Windows Server 2003 (Enterprise or Standard editions with latest service packs)
- Windows 7 Professional/Ultimate (32 or 64 Bit)
- Windows 8 Professional/Enterprise (32 or 64 Bit)
- Windows Server 2008 SP2 or higher (32 Bit)
- Windows Server 2008 R2 (64 Bit)
- Windows Server 2012 (64 Bit)